



EyeQue®



Remote
Control

What's in the Box



- 1 Back and Next buttons
- 2 Bluetooth indicator
- 3 Direction buttons
- 4 Confirm/ Select
- 5 Battery Indicator
- 6 On/Off Switch
- 7 Charging cable

Charging

- Charge the device using the micro USB cable provided (USB wall charger not included).

- LED on the battery indicator turns red while connected to a power source.

Safety

- This device contains a Lithium-ion battery.
- Do not discard in trash.
- Do not expose to high heat or incinerate.
- Do not charge in high temperatures.
- This device should be disposed of in accordance with your local rules and regulations.

To Use

- Turn on the Remote Control by using the on/off button.
- Bluetooth indicator light will blink blue indicating that the device is ready to be paired.
- Follow app instructions. The app will look for the Remote Control for up to 30 seconds.
DO NOT PAIR USING PHONE SETTINGS.

- Once the app discovers the Remote Control, it will pair automatically. The blinking blue light will remain solid when paired.

To Clean

- To clean, wipe the surface with a soft, dry cloth. Disinfect regularly with rubbing alcohol. Do not use soap or submerge in water.

EyeQue®
at-home vision tests

Made in China
Manufactured for
EyeQue Corp.
39608 Eureka Dr.
Newark CA
94560 USA

EyeQue.com

P/N 180-617-04
Rev 1.0